


Introduction

Gene therapy production is traditionally reliant on cells

AAV / viral vectors (encapsulated DNA)

- DNA payload and capsid delivery vector produced simultaneously intracellularly
- Technology Limitations
 - Long production times (weeks/months)
 - Low yields (10^{14} vg/L or <math><0.25</math> mg/L (DNA))
 - Propensity for product quality difficulties
 - Mix of partial, empty, full capsid vectors
 - Small-to-large scale variability
 - Production barriers (consistency, yields, and timelines) restrictive to stable preclinical / clinical development


Baculovirus / Sf9 facilitated AAV




Ohnori, T. "Advances in gene therapy for hemophilia: basis, current status, and future perspectives." *WJ* / *Hemophilia* 11(1), 31-41 (2020).

Non-viral (capsid free DNA)

- LNP mediated delivery
- Separates DNA drug substance from delivery vector for step wise control
 - Eliminates cell controlled encapsidation



GBIO v1 (old) DNA production – Baculovirus / Sf9 derived




How to address and overcome common gene therapy quality and manufacturing limitations?

GBIO v1 (old) non-viral DNA process


- Baculovirus / Sf9 facilitated production
- Standing complications:
 - Construct-to-construct variability
 - Extraneous DNA products
 - Low yields (1-10 mg/L)
 - Long production times (weeks/months)

Design

Rapid Enzymatic Synthesis (RES)




Learn more about iqDNA at Russell Mond's poster #1294




Experimental & Results

Characterizing and applying RES method



Demonstrating platform application



Objective & Approach

Develop a cell-free, enzymatic DNA production method

Goals

- Simplify system, remove complex (cellular & viral) components responsible for unrestrained byproduct formation [enhance process control]
- Eliminate cell growth times for rapid turnaround of DNA


Design Considerations

- Synthetic nature to facilitate heightened DNA functionality
 - Engineerable / modular design and components
- Platform
 - Build based on standard inputs and constituents
- Minimize number of enzymes to drive process simplicity
 - Fewer steps and components = robustness and reduced costs
- Empower scalability
 - Take advantage of highly processive nature of enzymes


Experimental & Results Continued

Accelerated and high-quality DNA production across scales

Process	Starting Scale (plasmid)			Amplification Scale (utilizes plasmid template from starting scale)					
	microgram - milligram			100 milligrams		1 gram			
	Plasmid receipt time	Processing time	Volumetric scale	Plasmid Amplification time	Processing time	Volumetric scale	Plasmid Amplification time	Processing time	Volumetric scale
RES	~2 wk	3 d	1-10 mL	2 d	4 d	~100 mL	2 d	4 d	~1 L
Sf9 Cell Based	~3 wk (bacmid)	3-4 wk	100-1000 mL	N/A	5-6 wk	>10 L	N/A	6-7 wk	>100 L



Expression validation using RES derived DNA



High quality, fully monomeric RES material displays enhanced expression compared with lower quality Sf9 material

Conclusion & Discussion

Rapid enzymatic synthesis (RES) demonstrated to be exceptional method for DNA gene therapy production

- Platform application
 - Standardized enzymes, recognition motifs, reaction conditions, and processing methods
- Synthetic makeup provides valuable engineerability and modular functionality
- Dramatically enhanced product quality control compared to cell-based methods
- Heightened R&D preclinical development enablement
 - At least 10-fold increase in screening / throughput compared with Sf9 process due to >100-fold reduction in volumetric scale
 - Plasmid receipt to purified drug substance in ~1 week (across all scales)
- Ultra-fast scalability at minimized production volumes
 - High volumetric (g/L) yields; traditional stirred tank reactor compatibility
 - ~1 week turn-around of gram quantities of drug substance with only $\mu\text{g}/\text{mg}$ starting plasmid required